



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



LID Hepta LED light (New)

- New,
- Made in France,
- A professional diagnostic light with a timeless design that combines power, robustness and versatility for intensive use.
- With its sturdy, flexible neck and ergonomic handle, it can be placed in any position and meets all medical lighting needs,
- An optional dimmer allows the operator to adjust the brightness for user comfort and optimal use,
- Specifications:
 - Light intensity: 980 Lumens,
 - Light color: 4500 K (neutral white),
 - Color rendering index: IRC 95/100,
 - Electricity consumption: 7 W,
 - Class 1 medical device,
 - Lifetime of LEDs minimum 50,000 hours of operation,
 - Electrical protection: insulation class II,
 - Standards: IEC 60 601 -1/ -2-41 IEC 60 598-2-25,
 - Power cord with a length of 3 m,
 - Available colors : white, black,
- Standard lamp length is 650 mm, for an extra charge you can purchase a lamp with a length of 1140 mm
- Dimensions :
 - Lamp head dimension: 100 x 110 mm,
 - Base dimension: 69 x 69 x 115 mm,
 - Weight: 1.6 / 2.4 kg,
 - Length: 650 / 1140 mm,
 - Interaxial spacing (4 holes in the lamp base): 54.5 mm,
- Diameter of the light field :
 - 30 cm distance: 130,000 lux field illumination 12 cm
 - Distance 50 cm: 60,000 lux field illumination 18 cm,
 - Distance 100 cm: 18 000 lux field illumination 30 cm,
- Additional mounting accessories are required for the lamp , of your choice (at an extra charge):
 - Telescopic stand base 60-100 cm metal white, heavy,
 - Telescopic stand base 60-100 cm white metal, heavy -version with brakes in 2 wheels,
 - Telescopic stand base 90-130 cm white metal, heavy,
 - Telescopic stand base 90-130 cm metal white, heavy -version with brakes in 2 wheels,
 - Plain clamp - black,
 - Clamp for medical rails with vertical lamp outlet,
 - Universal clamp,
 - Wall bracket,
 - ARM wall bracket white (2x40cm),
 - ARM wall mount white (1x40cm),
 - Tabletop mount ARM white (1x60cm),
 - Metal base white / gray,
 - Plate with metal pipe in the middle $\varnothing 16$,
 - Plate with metal tube on the edge $\varnothing 16$,
 - Plate on ARM holder,
 - Female socket with adaptor plate $\varnothing 16$ mm x 27 mm,
 - Adapter plate male rod $\varnothing 15$ mm x 30 mm,
- It has a valid Technical Passport issued,
- Warranty: 3 years,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

