



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



LID Carla LED light (New)

- New,
- Made in France,
- A treatment light with a modern and sophisticated design. The compact water drop-shaped head increases ease of use in many situations,
- Flexible and quiet boom provides excellent handling and precise and stable positioning,
- Very even, diffuse (25°) illumination,
- Specifications:
 - Light intensity: 420 Lumens,
 - Light color: 4500 K (neutral white),
 - Color rendering index: IRC 95/100,
 - Electricity consumption: 3 W,
 - Class 1 medical device,
 - Lifetime of LEDs minimum 50,000 hours of operation,
 - Electrical protection: insulation class II,
 - Standards: IEC 60 601 -1/ -2-41 IEC 60 598-2-25,
 - Power cord length: 3 meters,
 - Available colors : white, black,
- The standard length of the lamp is 650 mm, for an extra charge you can purchase a lamp with a length of 1000 mm or 1230 mm,
- Dimensions :
 - Light head dimension: 60 x 100 mm,
 - Base dimension: 69 x 69 x 60 mm,
 - Weight: 1.2 / 1.6 / 2.1 kg,
 - Length : 650 / 1000 / 1230 mm,
- Diameter of the light field :
 - 30 cm distance: 20,000 lux field illumination 15 cm
 - Distance 50 cm: 7,000 lux field illumination 25 cm,
 - Distance 100 cm: 1,700 lux field illumination 35 cm,
- Additional mounting accessories are required for the lamp , at your choice (at an extra charge):
 - Telescopic stand base 60-100 cm white metal, standard,
 - Telescopic stand base 60-100 cm metal white, standard -version with brakes in 2 wheels,
 - Telescopic stand base 90-130 cm metal white, standard,
 - Telescopic stand base 90-130 cm metal white, standard - version with brakes in 2 wheels,
 - Regular clamp - black,
 - Clamp for medical rails with vertical lamp outlet,
 - Universal clamp,
 - Wall bracket,
 - ARM wall bracket white (2x40cm),
 - ARM wall mount white (1x40cm),
 - Tabletop mount ARM white (1x60cm),
 - Metal base white / gray,
 - Plate with metal pipe in the middle Ø16,
 - Plate with metal tube on the edge Ø16,
 - Plate on ARM holder,
 - Female socket with adaptor plate Ø 16 mm x 27 mm,
 - Adapter plate male rod Ø 15 mm x 30 mm,
- It has a valid Technical Passport issued,
- Warranty: 3 years,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

