



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **LID Thera Infrared light (New)**

- New,
- Made in France,
- Infrared lamp mounted on a flexible gooseneck (length 86 cm, of which 34 cm is a rigid tube). Equipped with a handle, it can be easily set at the required distance,
- The model with a mechanical timer allows the operator to set the length of treatment from zero to 60 minutes with a beep when stopped,
- How to use: place the infrared lamp at a distance of about 40 cm: the degree of heat must always be pleasant. Sessions of 10 to 30 minutes are sufficient,
- Recommended for specialties such as physiotherapy / massage / acupuncture / orthopedics,
- Robust design allows for accurate and reliable positioning,
- Ergonomic handle allows precise placement,
- Specifications:
  - Electricity consumption: 250 W,
  - Wavelength: 1200 nm,
  - Class 1 medical device,
  - IR bulb life (E27): 5000 hours minimum,
  - Electrical protection: insulation class I,
  - Standards: IEC 60 601 -1/ -2-41 IEC 60 598-2-25,
  - Power cord length: 4 meters,
  - Available colors: white,
- Dimensions:
  - Dimension of the lamp head : 260 mm diameter,
  - Dimension of the base: 69 x 69 x 115 mm,
  - Weight : 2,9 kg,
  - Length : Flexible arm: 860 mm (including rigid 340 mm),
- The lamp can be retrofitted with a timer switch for an additional charge,
- Additional mounting accessories are required for the lamp , optional (for an extra charge):
  - Telescopic stand base 60-100 cm white metal, heavy,
  - Telescopic stand base 60-100 cm white metal, heavy -version with brakes in 2 wheels,
  - Telescopic stand base 90-130 cm white metal, heavy,
  - Stand telescopic base 90-130 cm metal white, heavy - version with brakes in 2 wheels,
  - Regular clamp - black,
  - Clamp for medical rails with vertical lamp outlet,
  - Universal clamp,
  - Wall bracket,
  - ARM wall bracket white (2x40cm),
  - ARM wall mount white (1x40cm),
  - Tabletop bracket ARM white (1x60cm),
  - Plate with metal pipe in the middle  $\varnothing 16$ ,
  - Plate with metal tube on the edge  $\varnothing 16$ ,
  - Plate on ARM holder,
  - Female socket with adaptor plate  $\varnothing 16$  mm x 27 mm,
  - Adapter plate male rod  $\varnothing 15$  mm x 30 mm,
- It has a valid Technical Passport issued,
- Warranty: 3 years,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

