



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Telescopic stand base - standard (New)**

- New,
- Dedicated for Tina and Carla lights,
- Made in France,
- Base in the shape of a five-pointed star with a diameter of 53 cm with wheels (diameter 5 cm). Very stable,
- Frame made of easy to clean and disinfect smooth steel coated with epoxy resin. Castered base without indentations or grooves, facilitating effective cleaning and reducing the retention of dust or other debris (ideal for hospitals and health centers),
- Adjustable column height, allowing the lamp to be adjusted to the appropriate height for all medical applications,
- Standard stand height is 60-100 cm , for an extra charge you can purchase a 90-130 cm stand,
- For an extra charge, it is possible to purchase a stand with additional brakes in 2 casters,
- Height: 60 - 100 cm / 90 - 130 cm ,
- Weight: 3.7 / 4 kg,
- Available colors : white or black,
- Warranty: 12 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

