



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Universal Clamp (New)**

- New,
- Made in France,
- Used for mounting on medical rails, tubes and table edges,
- Supplied with an adapter plate for the LID lamp, which can be placed on the clamp head or on the side of the clamp with 4-axis rotation (90°/180°/270°/360°),
- Allows the operator to position the lamp horizontally or vertically. Horizontal position is reserved only for Carla Focus L65cm and L100cm lamp,
- Maximum opening  $\varnothing$  4 cm,
- Made of anodized aluminum,
- With mounting plate to mount the lamp,
- Weight: 560 g,
- Warranty: 12 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

