



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Video Gastroscope PENTAX EG-2901 (Reconditioned)

- Condition: Reconditioned - Very Good,
- Made in Japan,
- Compatible with processor Pentax 3300,
- Direction of View: Forward,
- Outer diameter: 9,8 mm,
- Outer diameter insertion tube: 10,1 mm,
- Working field diameter: 2,8 mm,
- Working length: 1050 mm,
- Max Angulation Up/Down: 180°/90°
- Max Angulation Right/Left: 100°/100°,
- Field of View: 100°,
- Additional Equipment:
 - Original suitcase: 62 x 44 x 14 cm,
 - Valves,
 - Cover,
- Technical Passport,
- Warranty: 30-days

