



Grooming Table WESSELING-BV STABILO ORIGINAL (New)

- New,
- Made in the Netherlands,
- The frame is epoxy coated,
- Adjustable height: 38 – 97 cm,
- This couch is only available in electric height adjustable version, and one can choose a hand or footswitch (extra paid),
- This table has twin hinges and therefore it is very rigid,
- Max load: 200 kg,
- Dimensions:
 - Slip-resistant mat: 120 x 65 cm,
- **Additional Equipment (extra paid):**
 - an exhaust system with an opening,
 - handless for safety straps (random quantity and place),
 - connection of both controllers by splitter,
 - grooming handle (single or double),
 - castors (4 wheels with central lock),
 - sharp handle,
 - plastic container on a side,
 - stainless steel table top with mattress,
- Cat. No. 02,
- Technical Passport,
- 24-month warranty

