



Grooming Table WESSELING-BV STABILO COMPACT (New)

- New,
- Made in the Netherlands,
- Adjustable height: 45 – 100 cm,
- This couch is only available in electric height adjustable version, and one can choose a hand or footswitch (extra paid),
- Due to the scissor construction it has an upright height adjustment, and no sideways movement,
- Max load: 200 kg,
- Dimensions:
 - Stainless steel top: 105 x 55 cm,
 - Mattress: 100 x 50 cm,
 - Slip-resistant mat: 100 x 50 cm,
- **Additional Equipment (extra paid):**
 - an exhaust system with an opening,
 - handless for safety straps (random quantity and place),
 - connection of both controllers by splitter,
 - grooming handle (single or double),
 - castors (4 wheels with central lock),
 - sharp handle,
 - plastic container on a side,
- Cat. No. 03,
- Technical Passport,
- 24-month warranty



