



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



N.I.L. EVA sp1 (New)

Infrasonic Liposculpture

N.I.L. ® System EVA sp 1

EVA sp®1 is the basic system of the EVA sp® range. It regulates the compressed air supply to the Lipomatic® handpiece which allows the Nutational Infrasonic Liposculpture®.

System EVA SP 1:

- Lipomatic®3 or Lipomatic®4 handpiece,
- EVA sp®1 control box,
- Complete set of cannulas:
 - Mercedes 4 x 300,
 - 3G Accelerator 3 x 250,
 - 3G Accelerator 3 x 150,
 - Rebelo 3,5 x 250,
 - Infiltration 2,5 x 200,
 - Infiltration 2,5 x 300,
- Accessories:
 - 5x sewer jetter,
 - Regulator Key,
 - Cleaning Set,
 - Lipomatic ® Oil (125 ml),
 - Power supply unit,
 - Pedals.

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPAÑOL

WIR SPRECHEN DEUTSCH



