



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **REHABILITATION BED DREAM 321 (New)**

New  
Polish production  
The bed has a four-segment bed  
Autocontour function of the back and thigh segments  
Height adjustment in the range from 40 cm to 80 cm  
All functions of the bed are electrically controlled by means of the remote control  
Possibility of locking individual positions of the bed (Linak)  
Metal parts powder coated  
The surface of the lair is filled with a metal mesh  
Extra bed options (extra charge):  
Side barriers  
Boom with triangular handle  
Trendelenburg position  
Filling the fronts with a waterproof HPL board  
Four anti-static wheels with the possibility of blocking  
Technical parameters:  
LINAK drive and control  
Safe working load: 175 kg  
Maximum user weight: 135 kg  
Total weight: 92.6 kg  
External dimensions: length 216 cm, width 94 cm  
Length of the backrest segment: 80 cm  
Adjustment of the backrest segment: 0 ÷ 70 °  
Thigh rest segment adjustment: 0 ÷ 45 °  
Changing the height of the bed using electric motors: 40 - 80 cm  
Ground clearance under the bed - 22 cm  
Wheels 100 mm - standard with brakes  
The color of wooden board elements - beech  
The color of metal elements - RAL 7035  
Has a current Technical Passport issued  
Warranty 24 - months

