



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



REHABILITATION BED DREAM 321 (New)

New
Polish production
The bed has a four-segment bed
Autocontour function of the back and thigh segments
Height adjustment in the range from 40 cm to 80 cm
All functions of the bed are electrically controlled by means of the remote control
Possibility of locking individual positions of the bed (Linak)
Metal parts powder coated
The surface of the lair is filled with a metal mesh
Extra bed options (extra charge):
Side barriers
Boom with triangular handle
Trendelenburg position
Filling the fronts with a waterproof HPL board
Four anti-static wheels with the possibility of blocking
Technical parameters:
LINAK drive and control
Safe working load: 175 kg
Maximum user weight: 135 kg
Total weight: 92.6 kg
External dimensions: length 216 cm, width 94 cm
Length of the backrest segment: 80 cm
Adjustment of the backrest segment: 0 ÷ 70 °
Thigh rest segment adjustment: 0 ÷ 45 °
Changing the height of the bed using electric motors: 40 - 80 cm
Ground clearance under the bed - 22 cm
Wheels 100 mm - standard with brakes
The color of wooden board elements - beech
The color of metal elements - RAL 7035
Has a current Technical Passport issued
Warranty 24 - months

