



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Rehabilitation Bed DREAM 337 (Linak) (New)

- New
- Made in Poland
- 4 sections
- Auto Contour function for backrest / legplates
- Adjustable height: 40 cm to 80 cm
- Remote Control
- Position lock/block thanks to (Linak)
- Metal elements powder coated
- Bed bottom made from resilient panels
- **Additional Equipment (paid):**
 - Side rails
 - Trapeze stand
- Fronts fulfilment made from HPL board
- 4 casters with breaks
- Technical Data:
 - Power and control panel: Linak
 - Safe working load: 175 kg
 - Lifting capacity: 135 kg
 - Weight: 102,5 kg
 - Dimensions: L 215 cm, W 102 cm
 - Backrest length: 80 cm
 - Adjustable angle backrest: 0 ÷ 70°
 - Adjustable angle leg plates: 0 ÷ 45°
 - Electrical height adjustment: 40 - 80 cm
 - Clearance under the bed - 22 cm
 - Castors 100 mm
 - Wooden elements color - beech
 - Metal elements color according to - RAL 9006
- Technical Passport
- 24 months warranty

