



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Rehabilitation Bed DREAM 521 (New)

- New
- Made in Poland
- 4 sections
- 2 powerful actuators individually controlled
- Adjustable height: 31 - 71 cm
- Remote Control
- Metal elements powder coated
- Bed bottom made from resilient panels
- Technical Data:
 - Power and control panel: Linak
 - Safe working load: 175 kg
 - Lifting capacity: 135 kg
 - Weight: 78 kg
 - Dimensions: L 190 cm, W 90 cm
 - Backrest length: 80 cm
 - Adjustable angle backrest: 0 ÷ 70°
 - Adjustable angle leg plates: 0 ÷ 45°
 - Metal elements color according to - RAL 9006
- Technical Passport
- 24 months warranty

