



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Bedside Table DREAM 3 (New)

- New
- Made in Poland
- Made from particle board covered with veneer
- Wooden elements color: beech
- Frame made from powder coated steel
- On casters, two with central breaking/lock
- Technical Data:
 - Adjustable angle of countertop tilt (5 lvl): 0° - 45°
 - Adjustable height (6 lvl): 770mm - 1010 mm
 - External dimensions: 900 x 450 mm
 - Regulated countertop dimensions: 620 x 450 mm
 - Castors: 4/ 50 mm
 - Max lifting capacity: 20 kg
 - Weight: 18 kg
- Technical Passport
- 12 months warranty

