



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



BLOODY 02 blood collection chair (New)

New,

Made in Poland,

The design makes it easy to carry out the injection procedure on the left or right forearm of the person subject to the procedure, It is possible to order upholstery in a colour of your choice free of charge, except for "PRESTIGE" colours which are available at an extra charge,

The armchair is used in doctor's surgeries, health centres, hospitals and laboratories,

The armchair has a stable steel construction powder-coated, adjustable backrest to the lying position and two armrests adjustable in all planes,

The armchair can be fitted with

headrest,

adjustable footrest,

Technical parameters:

Height 100 cm,

Width 54 cm,

Depth 80 cm,

Seat height 50 cm,

Total length when unfolded 143 cm,

Weight 30 kg,

Maximum load 150 kg,

Headrest height adjustment 15 cm

Height adjustment of armrests 25 cm

Backrest angle adjustment 0°- 85°,

Headrest angle adjustment 0°- 90°,

It has a current Technical Passport,

Warranty: 24 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

