



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical Display BEACON C22S+ (New)

- New,
- 2MP Color Medical Display,
- Features:
 - MULTI,
 - BSS,
 - ALAS,
 - DICOM,
 - COLOR,
- Application:
 - CR/DR,
 - CT,
 - MRI,
 - PET / CT,
 - X-RAY,
 - DSA,
 - DSI,
 - 3D,
 - Nuclear Medicine,
 - CAD medical software,
 - Image reconstruction,
- Technical Specifications:
 - Panel Size: 21,3",
 - Support Color: 1,07 Billions,
 - Resolution: 1600 x 1200,
 - Pixel Pitch: 0,270 x 0,270 mm,
 - Brightness: 900 cd / m2,
 - Contrast: 1400 : 1,
 - Angle of View: 176°,
 - response Time: 20 ms,
 - Input:
 - DVI-D: 24 pin,
 - DisplayPort: 20 pin,
 - VGA: D-Sub 15 pin,
 - Dimensions: 369 (W) x 507 ÷ 627 (H) x 235 D) mm,
 - Weight: 9,2 kg,
 - Standard Mounting: VESA 100 x 100 mm,
 - Power supply: DC 12 V,
 - Max power consumption: 60 W,
 - Power Saving: EPA, VESA,
 - MTBF: > 50 000h,
 - Certificates: EN 60601, CCC, FCC, CE, UL,
- Technical Passport,
- 36-month warranty

