



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical Display BEACON G32S+ (New)

- New,
- 3MP Grayscale Medical Display,
- Features:
 - MULTI,
 - BSS,
 - ALAS,
 - DICOM,
 - QA,
- Application:
 - CR / DR,
 - CT,
 - MRI,
 - X-RAY,
 - DSA,
 - DSI,
- Technical Specifications:
 - Panel Size: 21,3",
 - Support Color: 4096 greyscale,
 - Resolution: 2048 x 1536,
 - Pixel Pitch: 0,2115 x 0,2115 mm,
 - Brightness: 1700 cd / m²,
 - Contrast: 1400 : 1,
 - Angel of View: 176°
 - Response Time: 20 ms,
 - Input:
 - DVI-D: 24 pin,
 - DisplayPort: 20 pin,
 - Dimensions: 369 (W) x 507 ÷ 627 (H) x 235 (D) mm,
 - Weight: 9,2 kg,
 - Standard Mounting: VESA 100 x 100 mm,
 - Power Supply: DC 12 V,
 - Max power consumption: 37 W,
 - Power saving: EPA, VESA,
 - MTBF: > 50 000h
 - Certificates: EN 60601, CCC, FCC, CE, UL
- Technical Passport,
- 36-month warranty

