



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Grayscale Medical Display BEACON G52S+ (New)

- New,
- 5MP Grayscale Medical Display,
- Features:
 - MULTI,
 - ALAS,
 - QA,
 - BSS,
 - DICOM,
- Application:
 - Digital Mammography,
 - DR,
- Technical Specifications:
 - Display Size: 21",
 - Support Color: 4096 grayscale,
 - Resolution: 2560 x 2048,
 - Pixel Pitch: 0,165 x 0,165 mm,
 - Brightness: 1200 cd/m2,
 - Contrast: 1200 : 1,
 - Angle of View: 170° / 170°
 - Response Time: 25 ms,
 - Input:
 - DVI-D: 24 pin,
 - DisplayPort: 20 pin,
 - Dimensions: 398 (W) x 510 ÷ 630 (H) x 235 (D) mm,
 - Weight: 10,5 kg,
 - Standard mounting: VESA 100 x 100 mm,
 - Power supply: DC 12 V,
 - Max power consumption: 47 W,
 - Power saving: EPA, VESA,
 - MTBF: > 50 000h,
 - Certificates: EN 60601, CCC, FCC, CE, UL,
- Technical Passport,
- 36-month warranty

