



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical Display BEACON G11S (New)

- New,
- 1MP Greyscale Medical Display,
- Features:
 - DISPLAY,
 - PERSPEKTIVE,
 - BSS,
 - DICOM,
- Application:
 - DSA,
 - DSI,
 - CT,
 - MRI,
- Technical Specifications:
 - Panel Size: 19",
 - Support Color: 1024 greyscale,
 - Resolution: 1280 x 1024,
 - Pixel Pitch: 0,294 x 0,294 mm,
 - Brightness: 1300 cd/m²,
 - Contrast: 1000 : 1,
 - Angle of View: 170°
 - Response Time: 30 ms,
 - Input:
 - DVI-D: 24 pin,
 - VGA: D-Sub 15 pin,
 - YPbPr: BNC*3,
 - Output: YPbPr: BNC*3
 - Dimensions: 417 (W) x 344 ÷ 410 (H) x 220 (D) mm,
 - Weight: 7,6 kg,
 - Standard Mounting: VESA 100 x 100 mm,
 - Power supply: AC 230 V,
 - Max power consumption 60 W,
 - Power saving: EPA, VESA,
 - MTBF: > 50 000 h,
 - Certificates: CCC, CB, FCC, CE, cTUVus,
- Technical Passport,
- 36-month warranty

