



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical Display BEACON E150 (New)

- New,
- 15" HD Endoscope & Surgical Display,
- Features:
 - MULTI,
 - COLOR,
 - DUAL,
 - IN/OUT,
 - BSS,
- Features:
 - Thoracic Endoscope / Laparoscopic Surgery,
 - Endoscope / Gastroscopy,
 - Operating Room,
 - ICU (Intensive Care Unit),
- Technical Specifications:
 - Panel Size: 15",
 - Support Color: 16,7 Billions,
 - Resolution: 1024 x 768,
 - Pixel Pitch: 0,297 x 0,297 mm,
 - Brightness: 350 cd/m²,
 - Contrast 600 : 1,
 - Angle of View: 170°,
 - Response Time: 8 ms,
 - Input:
 - DVI-D: 24 pin,
 - VGA: D-Sub 15 pin,
 - CVBS: BNC,
 - S-VIDEO: DIN4,
 - YPbPr: BNC*3,
 - RGBS: BNC*4,
 - Output:
 - DVI-D: 24 pin,
 - CVBS: BNC,
 - S-VIDEO: DIN4,
 - YPbPr: BNC*3,
 - RGBS: BNC*4,
 - Dimensions: 364 (W) x 274 - 364 (H) x 255 (D.) mm,
 - Weight: 7 kg,
 - Standard Mounting: VESA 100 x 100 mm,
 - Power supply: AC 230 V, / 50 - 60 Hz,
 - Max power consumption: 60 W,
 - Power saving: EPA, VESA,
 - MTBF: > 50 000 h,
 - Certificates: EN 60601, CCC, FCC, CE, UL,
- Technical Passport,
- 36-month warranty

