



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Medical Display BEACON E190 (New)**

- New,
- 19" HD Endoscope & Surgical Display,
- Features:
  - MULTI,
  - COLOR,
  - DUAL,
  - IN/OUT,
  - BSS,
- Application:
  - Thoracic Endoscope / Laparoscopic Surgery,
  - Enteroscope / Gastroscope,
  - Operating Room,
  - ICU (Intensive Care Unit),
- Technical Specifications:
  - Panel Size: 19",
  - Support Color: 16,7 Billions,
  - Resolution: 1280 x 1024,
  - Pixel Pitch: 0,294 x 0,294 mm,
  - Brightness: 300 cd/m<sup>2</sup>,
  - Contrast: 1000 : 1,
  - Angle of View: 178°,
  - Response Time: 25 ms,
  - Input:
    - DVI-D: 24 pin,
    - VGA: D-Sub 15 pin,
    - CVBS: BNC,
    - S-VIDEO: DIN4,
    - YPbPr: BNC\*3,
    - RGBS: BNC\*4,
  - Output:
    - DVI-D: 24 pin,
    - CVBS: BNC,
    - S-VIDEO: DIN4,
    - YPbPr: BNC\*3,
    - RGBS: BNC\*4,
  - Dimensions: 419 (W) x 441 - 564 (H) x 251 (D) mm,
  - Weight: 11 kg,
  - Standard Mounting: VESA 100 x 100 mm,
  - Power supply: DC 12 V,
  - Max power consumption: 45 W,
  - Power saving: EPA, VESA,
  - MTBF: > 50 000 h,
  - Certificates: EN 60601, CCC, FCC, CE, UL,
- Technical Passport,
- 36-month warranty

