



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical Display BEACON M652T (New)

- New,
- Color Medical Image Consultation Center,
- Features:
 - LARGE FORMAT,
 - BSS,
 - DICOM,
 - CONSISTENT,
- Application:
 - Comprehensive consultation,
 - Surgery teaching,
 - Training,
 - Telemedicine consultation,
- Dane techniczne:
 - Panel Size: 65 ",
 - Support Color: 1,07 Billions
 - Resolution: 3840 x 2160,
 - Pixel Pitch: 0,372 x 0,372 mm,
 - Brightness: 450 cd/m²,
 - Contrast 1400 : 1,
 - Angle of View: 178° vertical / 178° horizontal,
 - Response Time: 5 ms,
 - Input:
 - DisplayPort: 20 pin*4,
 - DVI-D: 24pin,
 - Dimensions: 1480 (W) x 890 (H) x 180 (D) mm,
 - Weight: 55 kg,
 - Power supply: AC 230 V,
 - Max power consumption: 200 W,
 - Power saving: EPA, VESA,
 - MTBF: > 50 000 h,
- Technical Passport,
- 36-month warranty

