



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Autoclave STATIM 2000 G4 (with Data Logger) (New)

- New,
- Made in Canada,
- Features:
 - Ultra fast 6:45 minute unwrapped cycle,
 - 10:40 minute wrapped cycle,
 - Exceptional instrument turnover,
 - 10 minute Dri-Tec drying Technology,
 - Single use water,
 - Aseptic transfer of instrument directly,
 - to the point-of-use,
 - Reduced instrument investment,
 - Gentle on instruments,
 - EN13060 compliant,
- Technical Data:
 - Class S,
 - Power input: 208 - 240 V,
 - Power: 1300 W,
 - Decibel rating - 56 dB,
 - Cycle registration- Data Logger,
 - Hollow wrapped 134 C - 10 mins,
 - Hollow unwrapped 134 C - 6 mins,
 - Rubber & Plastic 121 C - 20 mins,
 - Drying cycle minimum 60 mins,
 - Wrapped 72 mins,
 - Unwrapped - ca 36 mins,
 - A + B,
 - Cassette dimensions - 37 cm / 18 cm / 6,5 cm,
 - Cassette capacity - 1,8 l,
 - Reservoir capacity: 4 l,
 - Dimensions: 48 cm x 41 cm x 15 cm,
 - Menu in Polish,
 - Weight - 21 kg,
- Technical Passport,
- 24-month warranty

