



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Autoclave STATIM 5000 G4 WIFI (Data Logger built-in) (New)

New,
Canadian production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 1300 W,
Selectable menus in: Polish, English, Czech, Danish, Spanish, French, Hungarian, Italian, Korean, Dutch, Portuguese, Romanian, Russian, Ukrainian,
Advantages:
Colour touchscreen LCD display,
Advanced G4 operating system,
Easy to use,
Built-in Data Logger,
USB flash drive recording,
LAN network interface,
Energy saving,
Fast sterilization cycles - up to 9 min. for unpacked instruments,
Low consumption of distilled water,
Simple construction facilitating service and periodical maintenance,
Compliance with EN13060 for sterilisation of type A and B loads,
Technical data:
Class S,
Operating noise - 57 dB,
Cycle recording - built-in Data Logger,
Cycle for packed 134°C - 18 min,
Cycle time for unpacked 134°C - 9 min,
Plastics and rubber 121°C - 24 min,
Drying - minimum 60 min,
Sterilisation and drying of packaged - approx. 82 min,
Sterilisation and drying of non-packaged - approx. 45 min,
Sterilisation of batches of type - A and B,
Dimensions of the cassette (L/W/D) - 38 x 18 x 8 cm,
Capacity of the cassette - 5.1 l,
Water tank capacity - 4 l,
LAN operation : WiFi or network cable RJ45,
Equipped with RS232 port,
Pre-vacuum: thermodynamic,
Drying: thermal,
Included:
STATIM 5000-S cassette with shelf with drying accelerator plates,
STAT-DRI Plus 59ml with atomiser,
Test sample - Process Challenge Device,
Sterilisation emulsifiers - 20 pieces,
Waste bottle,
Outlet pipe,
Fittings for tube attachment,
USB memory stick 16GB,
USB memory stick with user manual in: Polish, Czech, Danish, Spanish, French, Hungarian, Italian, Korean, Dutch, Portuguese, Romanian, Russian, Ukrainian,
User manual in English,
Power cord,
Dimensions: 58 x 41 x 20 cm,
Set weight: 40 kg,
Holds a current Technical Passport,
Warranty: 12 months,
Finance options (Poland only): Instalments, Leasing, Loan

