



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



STATIMATIC + adapters (New)

- New,
- Made in Germany,
- Features:
 - Safe: validated cleaning verified by microbiological testing done by an independent laboratory.
 - Economical: low purchasing and operating cost,
 - Powerful: hard on germs, bacteria and dirt; gentle on your instruments,
 - User Friendly: very easy to use,
 - Space Saving: the STATIMATIC PLUS is designed to save valuable space in the hygiene room,
 - Environmentally Friendly: does not use damaging aerosol propellants,
 - Timing per cycle - 15 sec,
 - Auto detection,
 - Optional adapters: NSK, KaVo, Bien Air, W&H, Midwest,
 - CM Care protection system,
- Spray efficiency - 300 times,
- Capacity - 1-3,
- Noise level - 60 dB,
- Air pressure - 4 - 6 bar,
- Air demand - 40 NI / min,
- Power supply: 208 - 240 V,
- Dimensions: 19 cm x 29 cm x 40 cm,
- Weight - 8,6 kg,
- Technical Passport,
- 36-month warranty



