



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **STATIMATIC + adapters (New)**

- New,
- Made in Germany,
- Features:
  - The STATMATIC PLUS performs fully automated interior cleaning, which fully conforms with the strict German-based RKI guidelines\*. A water-based cleaner is used in order to facilitate the removal of bacteria, germs and protein. In accordance with the American medicinal cleaning norm, ASTM E2314-03, a germ/bacteria reduction of at least 2 log increments (i.e. >99% removal) is achieved. This process was tested using a microbiological method and was validated by mdt [Medical Device Testing GmbH], an independent, ZLG-accredited laboratory,
  - The process is divided into three phases:
    - 1. Cleaning of the water and air spray channels with STATMATIC PLUS cleaner,
    - 2. Cleaning of the drive air channel and hollow cavities with STATMATIC PLUS cleaner,
    - 3. Lubrication of all internal components with STATMATIC PLUS dental oil,
    - Each phase consists of several pulses during which the respective medium is applied to the inside of the handpiece. For contra-angles or slow-speed handpieces, a micro-motor rotates all components during the operation in order to ensure optimal distribution of the media,
    - Each phase is completed by a thorough purging with compressed air,
    - Timing - 15 sec,
  - Auto detection / CM Care,
  - Optional adapters: NSK, KaVo, Bien Air, W&H, Midwest,
  - Space Saving - the STATMATIC PLUS is designed to save valuable space in the hygiene room,
  - Environmentally Friendly - does not use damaging aerosol propellants,
  - User Friendly - very easy to use,
- Spray efficiency- 300 times,
- Tips qty - 1-3,
- Noise lvl - 60 dB,
- Air Pressure- 4 - 6 bar,
- Air supply - 40 NI / min,
- Power supply: 208 - 240 V,
- Dimensions: 19 cm x 29 cm x 40 cm,
- Weight - 8,6 kg,
- Technical Passport,
- 36-month warranty

