



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Sealing Device EUROSEAL 2001 Plus (New)

- New,
- Made in Italy,
- Features:
 - The sealing lever at the front makes packing easier,
 - The roll-locking bar keeps the roll stationary inside the sealing area at all times,
 - It is extremely ease to use and it limits the possibility of mistakes,
 - The high-quality materials and solid design make it indestructible,
 - The compact size allows it to be placed on a work surface or flush-mounted,
 - The smooth surfaces allow quick and thorough cleaning,
- Power supply - 230 V,
- Power - 110 W,
- Sleevs handle - 39 cm,
- Sealing width - 12 mm,
- Dimensions: 48 cm x 39 cm x 19 cm,
- Weight - 7,4 kg,
- Technical Passport,
- 12-month warranty

