



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Endoscopy table ALLOS M (New)

New,  
Made in Poland,  
Stable steel construction powder coated,  
Height adjustment of the table is carried out with a hand crank placed in the back part of the table frame,  
Standard features include two-part plywood sleep surface,  
Headrest is adjustable with a gas spring,  
There is a free option to order upholstery in the colour of your choice, except for "PRESTIGE" colours which are available at an extra charge,  
Table features and benefits  
Gas spring from the German company BANSBACH,  
Upholstery of very good quality,  
Foam, 5 cm thick,  
Steel ball bearings are used in the movable parts of the table, which makes the table very reliable,  
Frame elements made of the best Polish steel,  
As additional options are available (for an extra charge):  
Armrests,  
Hanger for pad,  
Castor system,  
Euro rails (40cm),  
Knee pads,  
Injection support,  
Technical parameters:  
Length 200cm  
Width 65/70/80 cm  
Height 52-90 cm  
Headrest angle adjustment -70° to 40°  
82 kg empty weight  
Maximum load 160 kg  
It has a valid Technical Passport,  
Warranty on structure 24 months,  
Warranty on hydraulic accessories 12 months,  
If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.  
WE SPEAK ENGLISH  
HABLAMOS ESPANOL  
WIR SPRECHEN DEUTSCH  
ON PARLE FRANCAIS  
WIJ SPREKEN NEDERLANDS  
МЫ ГОВОРИМ ПО РУССКИ  
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ  
WE SHIP WORLDWIDE

