



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Endoscopy Examination Table EOS E06 (New)**

New,  
Made in Poland,  
It has a stable steel construction coated with powder paint,  
Standard features include four-section sleep surface, foldable to Pivot and seat positions, made of plywood,  
Table height adjustment is carried out with an electric actuator operated by a manual remote control,  
Basic option - manual remote control - free of charge,  
Convertible option - foot remote control (no hand control) - for extra charge,  
Complete option - hand control + foot remote - at extra charge,  
Individual table top parts are adjustable with gas springs,  
There is a free option to order upholstery in the colour of your choice, except for "PRESTIGE" colours which are available for an extra charge,  
Table features and benefits  
Gas spring from German company BANSBACH,  
Electric actuator by a German company DEWERT,  
Upholstery of very good quality,  
5 cm thick foam,  
Steel ball bearings were used in the movable parts of the table, which makes the table very reliable,  
Frame elements made of the best Polish steel,  
As additional options are available (for an extra charge):  
Armrests,  
Hanger for the underlay,  
Castor system,  
Foot remote control (instead of manual remote control),  
Control from the frame,  
Euro rails (40cm),  
Knee bolsters,  
Injection support,  
Arm support for the chair,  
Technical parameters:  
Length 200cm,  
Width 65/70/80 cm,  
Height 55 - 106 cm,  
Trendelenburg position: 00 to 200,  
Anti-Trendelenburg position: - 180 to 00,  
Headrest angle adjustment: -75° to 45°,  
Adjustment of the footrest angle: - 450 to 500,  
Adjustment of the angle of the thoracic section: - 180 to 800,  
Lumbar section angle adjustment: -240 to 840,  
Weight 113 kg,  
Maximum load 160 kg,  
It has a valid Technical Passport,  
Warranty on structure 24 months,  
Warranty on hydraulic accessories 12 months,  
If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.  
WE SPEAK ENGLISH  
HABLAMOS ESPANOL  
WIR SPRECHEN DEUTSCH  
ON PARLE FRANCAIS  
WIJ SPREKEN NEDERLANDS  
МЫ ГОВОРИМ ПО РУССКИ  
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ  
WE SHIP WORLDWIDE

