



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Grooming Table WESSELING-BV STABILO SUPER (New)

- New,
- Made in the Netherlands,
- The frame is epoxy coated,
- Low maintenance hinges on the moving parts,
- All feet are adjustable for stability on uneven floors,
- This couch is only available in electric height adjustable version, and one can choose a hand or footswitch (extra paid),
- Technical Data:
 - Length: 120 cm,
 - Width: 65 cm,
 - Height Adjustment: 38 - 97 cm,
 - Max load: 175 kg,
- **Additional Equipment (extra paid):**
 - An exhaust system with an opening,
 - Handless for safety straps (random quantity and place),
 - Connection of both controllers by splitter,
 - Grooming handle (single or double),
 - Castors (4 wheels with central lock),
 - Sharp handle,
 - Plastic container on a side,
 - Stainless steel table top with mattress,
- Cat. No. 08,
- Technical Passport,
- 24-month warranty

