



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



BCI SpO2 Sensor (New)

- New,
- Made in Uruguay,
- Options:
 - Adult Clip,
 - Adult Silicone,
 - Pediatric Clip,
 - Pediatric Silicone,
 - Neonate/Multisite,
 - Wrap Silicone,
- **Cable Length 2,7 m** (cable extension available on demand),
- Compatible with:
 - BCI: 3300, 3301, 3302, 3303, 3304, 3401, 3403, 9200,
 - Spacelabs Burdick: PRO2,
 - Drager: Oxipac 2000,
 - MedChoice: MD300A,
 - Invacare 3300nv,
- **Cable Length 0,9m** (cable extension available on demand),
- Compatible with:
 - BCI: 3180, 3401, 3304, 3303, 3302, 3301, 3300, 3400, 3403,
 - Smiths Medical BCI Capnocheck 9004, Spectro 2 | 10, 20, 30, SurgiVet,
 - Spacelabs Burdick: PRO2,
 - Drager: Oxipac 2000,
 - Armstrong Medical AD-1000,
 - Beijing Choice Electronic: Pulse Oximeter model MD300I1.
- 12 months warranty,

