



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Criticare CSI Sp02 Sensor (New)

- New,
- Made in Uruguay,
- Options:
 - Adult Clip,
 - Adult Silicone,
 - Pediatric Clip,
 - Pediatric Silicone,
 - Neonate/Multisite,
 - Wrap Silicone,
- **Cable Length 2,7 m** (cable extension available on demand),
- Compatible with:
 - Scholar, Poet (602-3A), 1100, 2200, 500, 500A, 501, 502 series, 503, 504, 504/P, 504O/P, 504US/P, 507, 507E/P/R, 507/S, 507SD, 5070/P, 508, 509, 581,
 - Alaris 4415 (IVAC): VitalCheck 4400, 4410,
 - Criticare: 503DX-W, 503S, 503DX, 503 Spot, 504DX/P, 504DX, 506 series (except 506LN3), 507ELC2, 8000 series, 8100 series,
 - Poet Plus, Poet R Plus, Poet (602-3A),
 - Scholar, Scholar III,
 - nGenuity.
- **Cable Length 0,9m** (cable extension available on demand),
- Compatible with:
 - Alaris (IVAC): VitalCheck 4400 511SD,
 - Criticare: 503DX-W, 503S, 503DX, 503 Spot, Mini Spot, 504DX/P, 504DX, 506 series (except 506LN3), 507ELC2, 507EL, 8100 series, Poet Plus, Poet (602-3A),
 - Scholar, Scholar III.
- 12 months warranty,

