



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



GE-Datex Ohmeda SpO2 Sensor (New)

- New,
- Made in Uruguay,
- Options:
 - Adult Clip,
 - Adult Silicone,
 - Pediatric Clip,
 - Pediatric Silicone,
 - Neonate/Multisite,
 - Wrap Silicone,
- **Cable Length 2,7 m** (cable extension available on demand),
- Compatible with:
 - AS/3, AS/5, CS/3 Modular Monitors with Mseries modules, CH/S, CH/2S, CH/RS, CG/2S, Cardiocap I, Cardiocap II, Cardiocap 5, Cardiocap S/5 Light Monitor, Capnomac Ultima, M-Estpr, Oscar II, Satellite, Satellite Plus II, Trans Plus, ULTSVL, ULTSVO, ULTSIO, ULTS00, OSP-200.
 - S/5, i4 Modular Monitors con E-Series, Corometrics 259 with Trusignal Ohmeda® Technology.
 - Ohmeda Trusat,
 - Ohmeda 3700, 3710, 3770, 3775, 3800, 3900, Biox 3740, 3760, 3800, 4700 Oxicap, 5250 RGM, Sac II, Oxi-Tip Plus, Tuffsat, S/5.
- **Cable Length 0,9m** (cable extension available on demand),
- Compatible with:
 - Datex Ohmeda Trusignal,
 - Universal Oxytip,
 - AS/3, CS/3, CH/S, CH/2S, CH-RS, SAF-F, CG/ 2S, Cardiocap I, Cardiocap II, Cardiocap 5, Capnomac Ultima, Light Monitor, M-Estpr, Oscar II, OSP-200, Satellite Trans, Satellite, Satellite Plus II, Trans Plus, System 5, ULTS00, ULTSIO, ULTSVI, ULTSVO.
- 12 months warranty,

