



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



GE-Marquette SpO2 Sensor (New)

- New,
- Made in Uruguay,
- Options:
 - Adult Clip,
 - Adult Silicone,
 - Pediatric Clip,
 - Pediatric Silicone,
 - Neonate/Multisite,
 - Wrap Silicone,
- **Cable Length 2,7 m** (cable extension available on demand),
- Compatible with:
 - EAGLE, DASH, SOLAR Series, GE Responder 2000, Responder 3000, TRAM Modular x50SL Series Nellcor® Technology,
 - Dash 2000, 3000, 4000, 5000, i/9500 (Tram 451N5/851N5). Solar 7000, 8000 e Eagle 1000, 3000, Smk SMV, Corometrics 250 with Masimo technology,
 - Corometrics 250cx and Marquette monitors with Oximax® technology,
 - Hellige Eagle 1000, 3000, SMU610, 611, SMK210, SMV, Series with Nellcor® (standard) Technology,
- **Cable Length 0,9m** (cable extension available on demand),
- Compatible with:
 - GE - Marquette devices with Nellcor® technology,
- 12 months warranty,

