



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Invivo Research SpO2 Sensor (New)

- New,
- Made in Uruguay,
- Options:
 - Adult Clip,
 - Adult Silicone,
 - Pediatric Clip,
 - Pediatric Silicone,
 - Neonate/Multisite,
 - Wrap Silicone,
- **Cable Length 2,7 m** (cable extension available on demand),
- Compatible with:
 - 4500, 4500 Plus, Millenia, Sat Scout Omega 1445, Omni-Trak TVS,
- **Cable Length 0,9m** (cable extension available on demand),
- Compatible with:
 - All Invivo Research devices with 6pin DB9,
 - 4500 Plus 2 con Tecnología Nellcor®,
- 12 months warranty,

