



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MDE SpO2 Sensor (New)

- New,
- Made in Uruguay,
- Options:
 - Adult Clip,
 - Adult Silicone,
 - Pediatric Clip,
 - Pediatric Silicone,
 - Neonate/Multisite,
 - Wrap Silicone,
- **Cable Length 2,7 m** (cable extension available on demand),
- Compatible with:
 - Escort 1 / E 300, EsCodert 100 OPT 11A, EsCodert 300A, Escort II OPT30, Prism, 20300, 20401, 20403,
- 12 months warranty,



Model	Options	Price
MED	Adult Clip	80000 PLN netto
	Adult Silicone	80000 PLN netto
	Pediatric Clip	80000 PLN netto
	Pediatric Silicone	80000 PLN netto
	Neonate/Multisite	80000 PLN netto
	Wrap Silicone	80000 PLN netto