



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEK Sp02 Sensor (New)

- New,
- Made in Uruguay,
- Options:
 - Adult Clip,
 - Adult Silicone,
 - Pediatric Clip,
 - Pediatric Silicone,
 - Neonate/Multisite,
 - Wrap Silicone,
- **Cable Length 2,7 m** (cable extension available on demand),
- Compatible with:
 - MP-100, MP-110, MP-400, MP-500, MP-600, MP-1000,
 - MP110 Plus, MP111,
 - PAMO I, PAMO II, PAMO III, MP 570T, MP 700, MP 800C,
 - MP 1000NT, MP 1300,
- **Cable Length 0,9m** (cable extension available on demand),
- Compatible with:
 - MP110 Plus, MP111, PAMO I, PAMO II, PAMO III, MP 570T, MP 700, MP 800C, MP 1000NT, MP 1300
- 12 months warranty,

