



PRAISTON Sp. z o.o.
ul. Górska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Mindray SpO2 Sensor (New)

- New,
- Made in Uruguay,
- Options:
 - Adult Clip,
 - Adult Silicone,
 - Pediatric Clip,
 - Pediatric Silicone,
 - Neonate/Multisite,
 - Wrap Silicone,
- **Cable Length 2,7 m** (cable extension available on demand),
- Compatible with:
 - PM 600, PM 6000 (New), PM 7000, 6201, PM 8000, PM 9000, PM 9000 Express, MEC 100, MEC 200, MEC 1000, MEC 1200, 809B, VS 800; Mennen Enmove 1000, Mercury,
 - Kontron: 7138, 7840, 7845,
 - Medec: MEDEC-1000, M-1000, M-1100,
 - Creative: PC9000, UP6000,
 - Edan: M3, M3A, M3B, M8, M8A, M8B, M9, M9A, M50, M80, iM8, iM8A, iM8B, iM9, iM9A, iM50, iM80,
 - Goldway, Contec and Danatech, Penlon Intermed PM 8000, Vital Lite, Vital Line, Vital Wave,
 - PM 5000, PM 6000, PM 8000,
 - BeneView T5, T8 and iPM-8900,
- **Cable Length 0,9m** (cable extension available on demand),
- Compatible with:
 - All Mindray devices with 7pin DB9,
- 12 months warranty,

