



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Nellcor SpO2 Sensor (New)

- New,
- Made in Uruguay,
- Options:
 - Adult Clip,
 - Adult Silicone,
 - Pediatric Clip,
 - Pediatric Silicone,
 - Neonate/Multisite,
 - Wrap Silicone,
- **Cable Length 2,7 m** (cable extension available on demand),
- Compatible with:
 - N20, N180, N6000, NPB40, NPB75, NPB3900, NPB190, NPB195, NPB4000,
 - N-100, N-200, N-1000 and N-2000 Preamp cable; MDE Escort E100, Drager VitAlert 3200,
 - N-595 and N-600; CAS Med 511 Apnea Monitor,
 - Cardell Veterinary 9403, 9500HD,
 - Colin BP-S510 and T105NXM,
 - Digicare LW600,
 - Dinamap Carescape V100, ProCare, Pro Series, ProSeries V2, Pro 300, Pro 1000, Pro 1000 V3, E-NSAT Modules for i4 Modular Monitors with Oximax option,
 - Fukuda Denshi DS-7100,
 - GE Procare,
 - Welch Allyn Spot Vital Signs 42NTB-E1, 42NT0-E1 LXI,
 - Colin: BP88/306,
 - Critikon Dinamap ProCare, Pro Series, Pro 300, Pro 1000, M-NSAT Modules,
 - Datascope: Accutorr Plus, Duo, Passport 2, Passport XG, Spectrum, Trio,
 - Nellcor: NPB-290, NPB-295, N-390, N-395, N-3000, N-6000, NPB-290, NPB-295, N-390, N-395, N-3000, N-6000,
 - Nihon Kohden: BSM-2300, BSM-4100, BSM-5100, BSM- 4104A, BSM-4114A,
 - Ohmeda enGuard CM5; Physio Control Lifepak 12,
 - N 45, N 65, N 85, NPB 40, NPB 75, NBP-195,
 - Colin; GE Dinamap Compact Dinamap MPS,
 - All Welch Allyn (with Oximax® technology),
- **Cable Length 0,9m** (cable extension available on demand),
- Compatible with:
 - All devices with Nellcor® Technology,
 - All devices with Oximax® Technology,
- 12 months warranty,



