



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Nonin SpO2 Sensor (New)

- New,
- Made in Uruguay,
- Options:
 - Adult Clip,
 - Adult Silicone,
 - Pediatric Clip,
 - Pediatric Silicone,
 - Neonate/Multisite,
 - Wrap Silicone,
- **Cable Length 2,7 m** (cable extension available on demand),
- Compatible with:
 - 8500, 8550, 8600, 8700, 8800, 2500, 2500A, 9840, 9600, 9700, 2120, 7500, Respironics 902M Plus, Lifesense and PureSAT,
- **Cable Length 0,9m** (cable extension available on demand),
- Compatible with:
 - 8500, 8550, 8600, 8700, 8800, 2500, 2500A, 9840, 9600, 9700, 2120, 7500, Respironics 902M Plus, Lifesense and PureSAT,
- 12 months warranty,

