



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Spacelabs SpO2 Sensor (New)

- New,
- Made in Uruguay,
- Options:
 - Adult Clip,
 - Adult Silicone,
 - Pediatric Clip,
 - Pediatric Silicone,
 - Neonate/Multisite,
 - Wrap Silicone,
- **Cable Length 2,7 m** (cable extension available on demand),
- Compatible with:
 - 90351-D/G, 90369, 9046, 90466, 90467, 90489, 90651-07/08, 90419-01, 90651A-07,
 - 90303B, 90351-0/6, 90465, 90466, 90489, 90651A-08, 90651-A, IM77, 90308, Scout,
 - 1050, 1600, 1700, 90367, 90369, 90469, 90496, 90651-08, Ultraview, 91220, 91369, mCare 300, only with Nellcor® technology,
- **Cable Length 0,9m** (cable extension available on demand),
- Compatible with:
 - All devices with Spacelabs Nellcor® Technology,
- 12 months warranty

