



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



BCI Sp02 Intermediary Cable (New)

- New,
- Made in Uruguay,
- Sensor connection / device connection:
 - DB9 / DB9 (BIOOX.BC.IC2)
- **Cable Length 2,0 m** (cable extension available on demand),
- Compatible with:
 - 78100-B, 78600, Vital Line, Vital Wave,
- 12 months warranty

