



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Critikon Sp02 Intermediary Cable (New)

- New,
- Made in Uruguay,
- Sensor connection / device connection:
 - DB9 / REDEL 7 PIN (BIOOX.CK.IC2)
- **Cable Length 2,0 m** (cable extension available on demand),
- Compatible with:
 - Dinamap 8700, 8710, 8720 Pulse, 8721, 9710, 9700, 8280 and Oxyshuttle,
- 12 months warranty,

