



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Mindray SpO2 Intermediary Cable (New)**

- New,
- Made in Uruguay,
- Sensor connection / device connection:
  - DB9 / REDEL 6 PIN (BIOOX.MI.IC2),
  - DB9 / REDEL 12 PIN (BIOOX.MI2.IC2),
  - DB9 / REDEL 7 PIN (BIOOX.MI3.IC2),
- **Cable Length 2,0 m** (cable extension available on demand),
- Compatible with:
  - PM 600, PM 6000 (New), PM 7000, 6201, PM 8000, PM 9000, PM 9000 Express, MEC 100, MEC 200, MEC 1000, MEC 1200, 809B, VS 800; Mennen Enmove 1000, Mercury; Kontron: 7138, 7840, 7845; Medec: MEDEC-1000, M-1000, M-1100; Creative: PC9000, UP6000; Edan: M3, M3A, M3B, M8, M8A, M8B, M9, M9A, M50, M80, iM8, iM8A, iM8B, iM9, iM9A, iM50, iM80; Goldway, Contec and Danatech, Penlon Intermed PM 8000, Vital Lite, Vital Line, Vital Wave,
  - PM 5000, PM 6000, PM 8000,
  - BeneView T5, T8 and iPM-8900,
- 12 months warranty,

