



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Nellcor SpO2 Intermediary Cable (New)

- New,
- Made in Uruguay,
- Sensor connection / device connection:
 - DB9 / DB9 7 PIN (BIOOX.NC.IC2),
 - DB9 / REDEL 12 PIN (BIOOX.NCL.IC2),
 - DB9 / 3M 14 PIN (BIOOX.NCM.IC2),
 - DB9 / 3M 14 PIN (BIOOX.NCN.IC2),
 - DB9 / DB9 (BIOOX.NCO.IC2),
- **Cable Length 2,0 m** (cable extension available on demand),
- Compatible with:
 - N20, N180, N6000, NPB40, NPB75, NPB3900, NPB190, NPB195, NPB4000,
 - N-100, N-200, N-1000 and N-2000 Preamp cable; MDE Escort E100, Drager VitAlert 3200,
 - N-595 and N-600; CAS Med 511 Apnea Monitor; Cardell Veterinary 9403, 9500HD; Colin BP-S510 and T105NXM; Digicare LW600; Dinamap Carescape V100, ProCare, Pro Series, ProSeries V2, Pro 300, Pro 1000, Pro 1000 V3, E-NSAT Modules for i4 Modular Monitors with Oximax option; Fukuda Denshi DS-7100; GE Procure; Welch Allyn Spot Vital Signs 42NTB-E1, 42NT0-E1 LXi,
 - Colin: BP88/306; Critikon Dinamap ProCare, Pro Series, Pro 300, Pro 1000, M-NSAT Modules; Datascope: Accutorr Plus, Duo, Passport 2, Passport XG, Spectrum, Trio; Nellcor: NPB-290, NPB-295, N-390, N-395, N-3000, N-6000, NPB-290, NPB-295, N-390, N-395, N-3000, N-6000; Nihon Kohden: BSM-2300, BSM-4100, BSM-5100, BSM- 4104A, BSM-4114A; Ohmeda enGuard CM5; Physio Control Lifepak 12,
 - N 45, N 65, N 85, NPB 40, NPB 75, NBP-195; Colin; GE Dinamap Compact Dinamap MPS; All Welch Allyn . (with Oximax® technology),
- 12 months warranty,

