



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Temperature Probe - Dräger-Siemens (New)

- New,
- Made in Uruguay,
- Options:
  - **Dräger 8000 Series**(Cable length 1,5 m):
    - Adult Skin (BIOTEMP.2001PT),
    - Pediatric Skin (BIOTEMP.2001PTP),
  - **Dräger Caleo Series**(Cable length 1,5 m):
    - Adult Skin (BIOTEMP.2200PT),
    - Pediatric Skin (BIOTEMP.2200PTP),
  - **Infinity Vista, Delta, Gamma, Kappa**(Cable length 3,0 m):
    - Adult skin (BIOTEMP.1100PT),
    - Pediatric skin (BIOTEMP.1100PTP),
    - Tympanic (BIOTEMP.1100T),
    - Adult esophageal/rectal (BIOTEMP.1100RA),
    - Pediatric esophageal/rectal (BIOTEMP.1100RP),
- **Cable extension available on demand,**
- 12 months warranty,

