



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Temperature Probe - Fischer & Paykel (New)

- New,
- Made in Uruguay,
- Options:
 - **IW 700, IW 800 & IW 900 Series**(Cable length 1,5 m - 2.0 m):
 - Adult skin (BIOTEMP.900PT),
 - **MR 700**(Cable length 1,5 m - 2.0 m):
 - Air 2 (BIOTEMP.900AD),
- **Cable extension available on demand,**
- 12 months warranty,

