



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



Temperature Probe - Nihon Kohden (New)

- New,
- Made in Uruguay,
- Options:
 - Adult skin (BIOTEMP.400PT),
 - Pediatric skin (BIOTEMP.400PTP),
 - Timpanic (BIOTEMP.400T),
 - Adult esophageal/rectal (BIOTEMP.400RA),
 - Pediatric esophageal/rectal (BIOTEMP.400RP),
- **Cable Length 3,0 m** (cable extension available on demand),
- Compatible with:
 - All Nihon Kohden devices with ¼ plug in,
- 12 months warranty,



Nihon Kohden		Modelo de sonda	
Temperatura probe	Descripción	Cable	Conector
BIOTEMP.400PT	Para uso en piel adulta	3,0 m	1/4"
BIOTEMP.400PTP	Para uso en piel pediátrica	3,0 m	1/4"
BIOTEMP.400T	Para uso tímpanico	3,0 m	1/4"
BIOTEMP.400RA	Para uso esofágico/rectal adulto	3,0 m	1/4"
BIOTEMP.400RP	Para uso esofágico/rectal pediátrico	3,0 m	1/4"