



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



Temperature Probe - Welch Allyn (New)

- New,
- Made in Uruguay,
- Options:
 - Adult skin (BIOTEMP.400PT),
 - Pediatric skin (BIOTEMP.400PTP),
 - Tympanic (BIOTEMP.400T),
 - Adult esophageal/rectal (BIOTEMP.400RA),
 - Pediatric esophageal/rectal (BIOTEMP.400RP),
- **Cable Length 3,0 m** (cable extension available on demand),
- Compatible with:
 - Propaq 100/200 Series, Atlas Propaq Encore, Propaq CS, Ultra Smart Protocol,
- 12 months warranty,



Welch Allyn		
Temperature probe	Options	Code
Adult skin (BIOTEMP.400PT)	Adult skin (BIOTEMP.400PT)	Adult skin (BIOTEMP.400PT)
Pediatric skin (BIOTEMP.400PTP)	Pediatric skin (BIOTEMP.400PTP)	Pediatric skin (BIOTEMP.400PTP)
Tympanic (BIOTEMP.400T)	Tympanic (BIOTEMP.400T)	Tympanic (BIOTEMP.400T)
Adult esophageal/rectal (BIOTEMP.400RA)	Adult esophageal/rectal (BIOTEMP.400RA)	Adult esophageal/rectal (BIOTEMP.400RA)
Pediatric esophageal/rectal (BIOTEMP.400RP)	Pediatric esophageal/rectal (BIOTEMP.400RP)	Pediatric esophageal/rectal (BIOTEMP.400RP)