



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Temperature Probe - YSI 400 (New)

- New,
- Made in Uruguay,
- Options:
 - Adult skin (BIOTEMP.400PT),
 - Pediatric skin (BIOTEMP.400PTP),
 - Tympanic (BIOTEMP.400T),
 - Adult esophageal/rectal (BIOTEMP.400RA),
 - Pediatric esophageal/rectal (BIOTEMP.400RP),
 - Air 1 (BIOTEMP.400A1),
 - Air 2 (BIOTEMP.400A2),
- **Cable Length 3,0 m** (cable extension available on demand),
- Compatible with:
 - CASMed: 750; Criticare: 506, 507, 507E, 507ELC, 507S, 508, 1100, 2200, 8100 series; Datascope: Passport 2, Spectrum, Spectrum OR, Passport XG, Trio; Datex Ohmeda: S/3, CS/3, S/5, Cardiocap 1, Cardiocap 2, Cardiocap 5, Light Monitor, Satellite Plus; Dräger: Delta, Vista, Gamma, Kappa; Invivo: MDE Prism; Mindray: PM 7000, PM 8000, PM 9000 Express, PM 9000, VS-800; Nihon Kohden: BSM-2300 Life Scope I, BSM-4100 Life Scope P, BSM-5100 Life Scope A, BSM-9510 Life Scope M, BSM-9800 Life Scope S.; Welch Allyn: Atlas
- 12 months warranty,

